



Surgical Marker (RRS 18-11110)

The RRS markers provide an accurate way to mark the anatomical sites for correct blood pressure monitoring or surgical sites. RRS surgical markers are designed also for marking the patient's skin so that the surgeon can easily and accurately mark the incision site prior to surgery.

Advantages

- ◆ Provide clean and accurate lines on the patient's skin
- ◆ Available with flexible graduated ruler for more accurate marking

Specifications of The RRS Surgical Marker

- ◆ Violet color ink
- ◆ Non-toxic
- ◆ Not easily removed by water and Povidone-iodine (Betadine)
- ◆ Alcohol can be used to remove mark on the skin (complete removal of the marker after bathing).
- ◆ Sterilized by Gamma irradiation
- ◆ V-pack packaging

Note

- ◆ Do not use and resterilize the RRS Surgical marker if the V-pack packaging has been opened.
- ◆ This product is a single-use marker and is used for one patient.

Storage

.The marker should be stored at temperatures between 10-30 °C and 30-80% Relative Humidity

